





## Accessories for SPELECRAMAN

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery SPELECRAMAN

Qty.	Order no.	Description	
1 pcs	CABSTAT1	<b>μStat 8000(P)/4000(P) one node cable connector for conventional electrodes</b>  Connects individual nodes of the μStat 8000(P)/4000(P) to conventional electrodes	
1 pcs	RAMANPROBE	<b>Raman Probe</b>  Reflection probe designed to be used with a single excitation 785 nm wavelength (up to 500 mW). Suitable to work with DropSens Raman Cell for Screen-Printed Electrodes or with any conventional Raman set-up.	

## Optional accessories

Order no.	Description	
<b>DropView SPELEC</b>	<b>DropView SPELEC Software</b>	
<p>DropView SPELEC is a Spectroelectrochemical software that controls SPELEC instrument, offering a perfect synchronization of the optical and electrochemical measurements, as well as advanced tools for data treatment.</p>		
<b>RAMANCELL</b>	<b>Raman Cell for Screen-Printed Electrodes</b>	
<p>Black teflon reflection cell for performing Raman Spectroelectrochemistry with screen-printed electrodes in combination with ref. RAMANPROBE.</p>		