





Accessories for 2.917.0010

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.917.0010

Qty.	Order no.	Description	
1 pcs	1.917.0010	917 Coulometer	
		<p>The 917 Coulometer is a compact titrator for coulometric water determination according to Karl Fischer. "Reduced to the max" – operating module, stirrer and pump are integrated in the 917 Coulometer for rapid reagent replacement. This saves a great deal of space on the laboratory table.</p> <p>Recommended measuring range: 10 µg...200 mg absolute water Coulometry is the ideal method for the determination of trace water in liquids, solids and gases. In addition, coulometry is an absolute method and thus no titer determination is necessary.</p>	
1 pcs	6.0341.100	Double Pt-wire electrode for coulometry	
		Indicator electrode used for coulometric Karl Fischer titration.	
1 pcs	6.0344.100	Generator electrode with diaphragm	
		<p>Generator electrode with diaphragm for coulometric Karl Fischer titrations, SGJ 29/22. This generator electrode is well suited for:</p> <ul style="list-style-type: none"> • samples with very low absolute water contents < 50 ppm • coulometric water content determination in oil samples • applications where ketone reagents are used • applications where low-conducting reagents are used (addition of chloroform, xylene >10%) 	
1 pcs	6.1403.030	KF adsorber tube for coulometer cell	
		For KF coulometers.	

1 pcs **6.1437.000** **SGJ stopper SGJ 14**



1 pcs **6.1439.010** **KF inlet and aspiration tube**
KF inlet and aspiration tube for coulometer cell



2 pcs **6.1446.040** **Threaded stopper / M6**



1 pcs **6.1446.090** **Stopper / B-14/15 / M8**



2 pcs **6.1448.020** **Septum 16 mm, 5 pieces**
Set of 5 pcs.



1 pcs **6.1464.320** **KF titration vessel / 80-250 mL / coulometric**



2 pcs **6.1602.105** **Bottle attachment / GL 45 / green**

For direct dosing from reagent bottles with GL 45 thread.



1 pcs **6.1608.023** **Amber glass bottle / 1000 mL / GL 45**

For Exchange Units. Bottle for auxiliary solutions.



1 pcs **6.1608.030** **Round glass bottle / 1000 mL / GL 45**



1 pcs **6.1609.010** **Adsorber tube with tubing nipple**

Lid with tubing nipple.



1 pcs **6.1623.000** **Overflow protection for Ti Stand**
For 803 Ti Stand



1 pcs **6.1801.120** **PVC tubing / 4 mm / 6 mm / 2 m**



2 pcs **6.1805.200** **PTFE tubing / M8 / 0.5 m**



1 pcs **6.1808.050** **Tubing adapter olive / M8 outer**
1 M8 outer thread and 1 tubing olive. E.g. for thermostat jacket of exchange units and stability measuring instruments.



1 pcs **6.1819.030** **Tube**
For M8 threaded holes. For 6.1602.1XX bottle attachments and 6.1608.0XX bottles



1 pcs **6.1819.050** **Tube**
For M8 threaded holes of 6.1602.10X Bottle attachments with 6.1608.030 Aspiration bottle. For 703 and 803 Ti Stand, 681 Pump Unit and 679 Rancimat



1 pcs **6.1903.030** **Stirring bar / 25 mm**
Stirring bar with magnetic core, PTFE covering, length 25 mm.



1 pcs **6.2023.020** **Clip for SGJ 14/15**
Clip for SGJ 14/16



1 pcs **6.2047.040** **Titration vessel holder for 899/917**
Titration vessel holder for coulometric cells for 899 Coulometer and 917 Ti-Touch



1 pcs **6.2104.020** **Electrode cable / 1 m / F**
For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).



1 pcs	6.2104.120	Electrode cable / 1 m / H	With plug H. For connecting generator electrode - KF Coulometer.	
1 pcs	6.2164.010	Power supply unit 100 - 240 V/24 V DC	Power supply unit for 915 KF Ti-Touch and 916 Ti-Touch	
1 pcs	6.2701.040	GL 18 screw cap with hole	Made of PBT. With a hole for the 6.1448.020 septum. For the 6.5405.000 cell without diaphragm and the 6.1464.32X and 6.1465.320 titration vessels.	
3 pcs	6.2713.000	Sleeve	For the fat-free ground-joint connection with the standard ground-joint SGJ 14.	
1 pcs	6.2713.010	Sleeve	Sleeve for the fat-free ground-joint connection with the standard ground-joint SGJ 29.	
1 pcs	6.2713.020	Sleeve	For the fat-free ground-joint connection with the standard ground-joint SGJ 19.	

1 pcs **6.2738.000** **Funnel**

For the KF coulometer.



1 pcs **6.2739.000** **Wrench**

For tightening connectors



1 pcs **6.2811.000** **Molecular sieve**

Molecular sieve. Bottle containing 250 g. Pore size: 0.3 nm. Without moisture indicator. For Rancimats and Karl Fischer instruments.



2 pcs **6.2816.030** **Needle with Luer connector**

Needle for filling in mercury into the Multi-Mode Electrode and for sample addition in Karl Fischer titration.



2 pcs **6.2816.090** **Syringe 5 mL**

5 mL syringe made of PP, with Luer connector.



1 pcs **6.6063.901** **Demo CD for tiBase 1.1**

Demo CD for tiBase. Test version for 30 days.



1 pcs **6.6066.000** **USB flash drive for Ti-Touch**

USB flash drive for 915 KF Ti-Touch and 916 Ti-Touch.



Optional accessories

Order no.	Description
-----------	-------------

2.142.0100 **Custom Q3X thermal printer**

Compact printer with USB interface for

- 900 Touch Control
- 915 KF Ti-Touch
- 916 Ti-Touch
- 917 Coulometer
- 877 / 848 Titrino plus
- 865 / 876 Dosimat plus
- 91X Meter (cable 6.2151.140)
- Eco Dosimat / Titrator
- 862 Compact Titrosampler
- 870 KF Titrino plus
- 899 Coulometer



Paper width 60 mm (40 characters). Including 6.2151.120 USB cable.

6.2723.300 **Protective cover for 900 Touch Control**

Protective cover for 900 Touch Control

