






## Accessories for 2.858.0110

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 2.858.0110

Qty.	Order no.	Description	
1 pcs	1.858.0020	<b>858 Professional Sample Processor – Pump</b>	
		The 858 Professional Sample Processor – Pump processes samples from 500 µL to 500 mL. The sample transfer takes place either with the installed bidirectional two-channel peristaltic pump or with an 800 Dosino.	
1 pcs	6.1803.020	<b>PTFE capillary 0.97 mm i.d. / 5 m</b> For all IC instruments	
1 pcs	6.1803.160	<b>PTFE capillary 0.75 mm i.d. / 3 m</b>  PEEK capillary for use with adsorber module.	
1 pcs	6.1815.010	<b>Spiral band 0.5 m</b> For holding together different cables or tubing.	
2 pcs	6.1826.020	<b>PVC pump tubing (blue/blue)</b> For online IC instruments and automation in voltammetry.	

---

**1 pcs**     **6.1831.020**     **PEEK capillary 0.75 mm i.d., 22 cm**  
For all IC components.



---

**1 pcs**     **6.1835.050**     **PEEK sample tube**  
For sample changers for transferring larger-sized sample volumes, e.g. for brightener determination.



---

**1 pcs**     **6.2041.450**     **Sample rack 56 x 11 mL + 56 x 50 mL**  
Multi-row sample rack for 56 samples in 11 mL sample vessels (6.2743.050) and 56 samples in sample vessels of 50 mL (6.2747.010).



---

**1 pcs**     **6.2151.000**     **Cable USB A – mini-DIN 8-pin**  
Controller cable



---

**1 pcs**     **6.2621.080**     **Capillary cutter**  
For plastic capillaries. Used with IC instruments



---

**1 pcs**     **6.2621.090**     **1/2 in. wrench**  
For 1/2 in. screws. For IC Sample Processors



---

**1 pcs**     **6.2621.100**     **3 mm hex key for IC Sample Processors**



---

**1 pcs**     **6.2621.130**     **Hexagon key 2 mm**  
2 mm.



---

**1 pcs**     **6.2739.000**     **Wrench**  
For tightening connectors



---

**1 pcs**     **6.2743.057**     **Sample tubes 11 mL**  
200 pieces. For IC Sample Processors and VA Autosampler



---

**2 pcs**     **6.2744.010**     **Pressure screw 5x**  
With UNF 10/32 connection. For the connection of PEEK capillaries



---

**2 pcs**     **6.2744.060**     **Threaded stopper**  
For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



---

**1 pcs**     **6.2744.080**     **M6 thread / UNF 10/32 adapter**  
With M6 outer thread. For IC instruments.



---

**2 pcs**     **6.2744.160**     **Nozzle/UNF 10/32 coupling with security device**  
Security device for connecting pump tubing to a nozzle.



---

**1 pcs**     **6.2744.290**     **T connector 3 × UNF 10/32 with holder**



---

**4 pcs**      **6.2747.010**      **Sample vessel 50 mL**

With blue screw cap and mL graduation. Set of 25 pieces.



---

**1 pcs**      **6.2751.110**      **Safety shield**

Safety shield for the Advanced IC/VA Sample Processors.



---

**1 pcs**      **6.2846.000**      **Sample needle made of zirconium oxide**

Sample needle for Sample Processors and Autosampler in ion chromatography and voltammetry.

