


## Accessories for 6.2041.340

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

## Optional accessories

Order no.	Description	
<b>6.1432.210</b>	<b>Sample beaker, clear glass, 75 mL, 1 piece</b>	
	Sample beaker, clear glass, 75 mL, suitable for sample racks: 6.2041.340, 6.2041.350, 6.2041.800 and 6.2041.900, 1 piece	
<b>6.1458.020</b>	<b>Titration head / 4 x M10 connectors</b>	
	With 4 M10 screw nipples, including stoppers.	
<b>6.1458.030</b>	<b>KF titration head / 4 x M10 connectors</b>	
	With 4 M10 screw nipples, including stoppers.	
<b>6.1458.070</b>	<b>Titration head / 2 x SGJ 14</b>	
	Titration head for 75 mL beakers. With 2 openings SGJ 14, 2 openings 6.4 mm, and 3 pre-mounted dosing tips. Suitable for rinsing and aspiration equipment.	
<b>6.1459.400</b>	<b>Sample beaker, plastic (PE), 75 mL, 144 pieces</b>	
	Sample beaker, plastic (polyethylene, PE) 75 mL, suitable for sample racks: 6.2041.340, 6.2041.350, and 6.2041.800, suitable for 6.02041.040 OMNIS Software sample rack, 144 pieces.	