



Accessories for 6.2041.470

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Optional accessories

Order no.	Description	
6.1459.300	Sample beaker, plastic (PP), 120 mL, 250 pieces	
	Sample beaker, plastic (polypropylene, PP), 120 mL, suitable for sample racks: 6.2041.470, 6.2041.840, 6.2041.850, 6.2068.010, 6.2068.020, 6.2068.030, and 6.9920.141, 250 pieces	
6.1909.050	Propeller for 120 mL beaker (6.1459.300)	
	Propeller for 802 Stirrer when this is used with sample racks for 120 mL beakers (6.1459.300)	