



## Accessories for 6.2816.020

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

## Optional accessories

Order no.	Description	
<b>6.2744.120</b>	<b>Luer F pressure screw adapter</b> Connection between pressure screw and syringe.	
<b>6.2816.030</b>	<b>Needle with Luer connector</b> Needle for filling in mercury into the Multi-Mode Electrode and for sample addition in Karl Fischer titration.	
<b>6.2816.040</b>	<b>Purge tubing</b> With PTFE tubing and Luer connection. For syringes. For aspirating eluents.	